




Oxford Expression Technologies Ltd
 Oxford Brookes University
 Gypsy Lane Campus
 Oxford OX3 0BP UK
 t: +44 (0)1865 483236
 f: +44 (0)1865 483250
 e: info@oetltd.com

pOET2 transfer plasmid

pOET2 is a baculovirus transfer vector designed for high level expression of foreign genes under the powerful AcMNPV polyhedron (*polh*) promoter. The vector is smaller than other available transfer vectors (4547bp) which greatly facilitate the cloning steps. It has a bacterial origin of replication and an ampicillin resistance gene for selection in *E. coli*. The *polh* sequences have been replaced by a multiple cloning site in the reverse orientation to pOET1, containing unique restriction sites for insertion of the foreign gene. The AcMNPV sequences flanking the gene in the transfer vectors MCS allow recombination with the viral DNA to insert the expression cassette into the *polh* locus. pOET2 is compatible with any baculovirus system that utilizes homologous recombination in insect cells

AcMNPV 623-1769
 Pr_PH 1770-1870
 MCS 1880-1962
 AcMNPV 1995-2579
 Col_E1 2785-3404
 Ampicillin 3559-4419



Multiple Cloning Site - 1880 to 1962

	NotI	XhoI								
	EagI	Asp718I		SmaI						
	BseX3I	KpnI		EcoRI	XmaI	SalI	HindIII			
	BglII	SacII	SacI	PstI	ApaI	XbaI	BamHI			
1880	AGATCTGCGGCCGCGGTACCTCGAGAGCTCGAATTCTGCAGCCCGGGGCCCCGTGACTCTAGAAGCTTGGATCC									
	Q I C G R G T S R A R I L Q P G G P S T L E A W I									
	PacI									
1955	TTAATTAA									
	L N Z									