


**OXFORD**  
EXPRESSION  
TECHNOLOGIES

# QUICK START GUIDE to pOET2.1 N/C\_6xHis

Catalogue Number	2001031
Storage	Tightly capped at -20°C
Product Guarantee	1 Year from the date of purchase, when properly stored and handled

pOET2N/C 6xHis is a baculovirus transfer vector designed for high level expression of foreign genes under the powerful AcMNPV polyhedrin (polh) promoter (Pr\_PH). The vector encodes an N-terminal 6xHis-Tag® fusion sequence that may be utilised if the insert includes a stop codon. This greatly eases the purification of the recombinant protein since the 6xHis-containing fusion proteins bind with high affinity to Ni-NTA Agarose. If required, the 6xHis-Tag® can be removed by incubating the fusion protein in the presence of the proteinase cleavage enzyme Thrombin. There is also a 6xHis-Tag® for C-terminal fusions where the insert's own start codon can be used to replace the start codon supplied in pOET2N/C. It has a Col E1 origin of replication and an ampicillin resistance gene for selection in E. coli. The polh sequences have been replaced by a multiple cloning site (MCS) containing unique restriction sites for insertion of the foreign gene in the correct orientation. The AcMNPV sequences flanking the gene in the transfer vector's MCS allow recombination with the viral DNA to insert the expression cassette into the polh locus. pOET2N/C is compatible with any baculovirus system that utilises homologous recombination in insect cells.



## Multiple Cloning Site

CATCATCACCACCATCACACCGGTCTGGTCCGCGTGGATCA

6xHis

thrombin

BglII NotI Asp718I XhoI SacI EcoRI PstI SmaI ApaI SalI XbaI HindIII BamHI  
EagI SacI KpnI

AGATCTGCGGCCGCGGTACCTCGAGAGCTCGAATTCTGCAGCCCGGGCCGTCGACTCTAGAAGCTTGGATCC

SpeI PaeI

CATCATCACCACCATCACTAGTTAATTA

6xHis

## Oxford Expression Technologies Ltd

Bioinnovation Hub, Gipsy Lane, Oxford. OX3 0BP

t: +44 (0) 1865 483236

f: +44 (0) 1865 483250

e: info@oetltd.com

w: www.oetltd.com

